METTLER TOLEDO

InPro 325X (i) pH Electrodes

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	4
de	Betriebsanleitung	6
en	Instruction Manual	8
es	Manual de instrucciones	10
fi	Käyttöohje	12
fr	Instructions d'utilisation	14
hu	Használati utasítás	16
it	Istruzioni per l'uso	18
ja	取扱説明書	20
ko	사용 매뉴얼	22
nl	Gebruikershandleiding	24
pl	Instrukcja obsługi	26
pt	Manual de instruções	28
ru	Инструкция по эксплуатации	30
sv	Bruksanvisning	32
th	คู่มือคำแนะนำ	34
zh	说明书	36

InPro 325X(i) pH elektrody

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1 Úvod

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Instruction manual

1 Introduction

The METTLER TOLEDO InPro™ 325X(i) pH single rod measuring cell is a low maintenance, autoclavable and pressure resistant pH electrode with a pre-pressurised liquid electrolyte and a reference electrolyte together with a silver ions barrier. Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use. We recommend that you only operate the electrode in conjunction with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

2 Safety instructions

When cleaning or calibrating the electrode using solutions containing acids or alkalis both protective spectacles and protective gloves should be worn. Do not touch the electrode at the plug-in contacts, since electrostatic discharges can damage the electronics.

3 Product description

The inscription on each electrode contains the following information:

METTLER TOLEDO	Manufacturer of the electrode
InPro 325X i/SG/yyyy/zzz	Type designation: X = membrane glass i = ISM Digital/ SG = auxiliary electrode/ yyy = rod length in mm/ zzz = type of temperature sensor
combination pH	Type of electrode
pH 0 ... 14 (InPro 3250)	pH-measurement range
pH 0 ... 12 (InPro 3253)	
pH 1 ... 11 (InPro 3251, InPro 3252)	
<ul style="list-style-type: none"> • 0 ... 100 °C (up to 140 °C InPro 3250, InPro 3253) • -25 ... 80 °C (InPro 3251) • 0 ... 80 °C (InPro 3252) 	Temperature range for operation (for sterilisation)
Order No. 52 00X XXX	Order number
 SEV14 ATEX 0168 X	Identification and number of certificate
 IECEx SEV 14.0025X	Identification and number of certificate
 CML 22 UKEX 2108X	Identification and number of certificate
IS/I, II, III/1/ABCDEFG/T6	FM tested
See instruction manual	Note

In addition each electrode is provided with a serial number on the red part of the plug-in head in order to enable identification.

This product complies with the ATEX guideline 94/9/EC. The complete declaration can be downloaded under www.mt.com/pro.

4 Installation and commissioning

1. While unpacking check the electrode for mechanical damage. Please report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the protective cap and rinse the electrode for a short time with deionised water. Carefully remove the silicon protection from the diaphragm with the aid of the knife supplied (cf. green instruction images on the packaging). If the electrode has already been used the sensitive part can be cleaned if necessary with a soft toothbrush and a mild soap solution, and then rinsed for a short time with deionised water. After rinsing the electrode should only be dabbed. Any rubbing of the pH-sensitive glass can increase the response time as a result of electrostatic charge accumulation.

3. Check the space behind the pH sensitive glass for air bubbles and remove any air bubbles that may be present by shaking the electrode lightly in a vertical direction.
4. Place the electrode into the armature as described in the instructions for the armature.
5. Connect the electrode and the pH transmitter with the corresponding cable (VP6, VP8, AK9). Observe the connections scheme that accompanies the cable or the transmitter.
6. Intelligent ISM™ sensors such as InPro 325Xi enable «plug and measure» and extended diagnostics. Operation of the electrodes in the ISM mode requires the M700 transmitter as well as an ISM pH-module (2700i X). For installation, commissioning and operation of the ISM system please observe the instructions for the transmitter, module and cable.

5 Operation:

Calibration of the electrode and pH transmitter

Before a calibration immerse the electrode for 10 minutes in a pH 7.00 or pH 4.01 buffer solution whilst connecting it to the transmitter. A 2-point calibration is recommended, e.g. pH 7.00 and pH 4.01 buffer solutions. For further details please consult the operating instructions for the pH transmitter.

6 Maintenance

1. After each working cycle the electrode tip and diaphragm(s) are to be carefully rinsed with deionised water. Any drying out of residues from the solution being measured on these parts must be avoided at all costs!
2. If the electrode is not in use it is stored with the electrode tip and the diaphragm immersed in a tank of electrolyte (9823/3M KCl or 9848 FriscoIy™).
3. If the electrode is stored in its armature the rules described under Item 2 apply, although the electrolyte has to be slightly modified for storage purposes (2 parts of buffer solution pH 9.2 to 10 parts of electrolyte) in order to avoid any corrosion of the armature parts.
4. If an electrode is inadvertently stored in a dry condition for several days it must be rinsed for several hours in a standard storage electrolyte.
5. The plug-in head should occasionally be checked for possible traces of moisture. If necessary clean it thoroughly with deionised water or alcohol and then dry it off carefully.

7 Removal of sources of malfunction

With a soft toothbrush and water any residues can easily be removed from the electrode. Mild washing agents can also be used. Protein contaminations on the diaphragm can be most likely removed with the cleaning solution (order no. 51 340 070).

 **Warning:** Solutions containing hydrofluoric acid damage the pH sensitive part of the electrode.

8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

9 Warranty

On manufacturing defects 12 months after delivery.

10 Storage conditions

We recommend that the electrodes are stored at room temperature.

InPro, ISM and FriscoIy are trademarks of the METTLER TOLEDO Group.

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本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：
Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

* For digital pH sensor only (InPro3XXXi)

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 08/2022 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 52 002 631 C

Management System
certified according to
ISO 9001 / ISO 14001

**UK
CA**₂₅₀₃

CE₁₂₅₈

EAC

www.mt.com/pro

* 5 2 0 0 2 6 3 1 C *